

# RESMED




#### PRODUCT CODE

##### ResMed Hospital Full Face Mask (vented)

Small	60700
Medium	60701
Large	60702
Headgear (10 per pack)	60710

##### ResMed Hospital Nasal Mask (vented)

Medium	61104
Large	61105
Headgear (10 per pack)	61118



# RESMED

## Disposable Masks

Quick to fit. Effective to use.




ResMed, a leading global supplier of masks, has developed a line of custom-designed products for hospital use.

**Quick to fit and easy to use.** ResMed's disposable masks and circuits provide **cost-effective, time-saving solutions** for the clinical environment. Ideally suited for CPAP, bilevel pressure treatment and non-invasive ventilation therapy. **Single patient use only. Maximum 7 days.**

ResMed Corp Poway, CA, USA +1 858 746 2400 OR 1 800 424 0737 (toll free), 1 Elizabeth MacArthur Drive Bella Vista NSW 2153 Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). Hospital Full Face Mask protected by patents: AU 744593, AU 775051, CA 2261790, EP 0956069, US 6112746, US 6374826, US 6513526, US 6581602, US 6634358, US 6701927. Other patents and design registrations pending. Hospital Nasal Mask patents and design registrations pending. Mirage is a trademark of ResMed Ltd. © 2006 ResMed Ltd. Specifications may change without notice. 1010610/2 06 06

## ResMed Hospital Full Face Mask


The proven comfort and seal of ResMed's Mirage™ cushion technology in a hospital mask. Quick to fit and easy to adjust, the full face mask is available in three sizes—**small, medium and large**—providing a comfortable fit for a wide range of patients. Designed for single patient use, the masks can be easily discarded, eliminating cross-contamination concerns and costs associated with mask disinfection and sterilization.



ResMed flow generator setting: Mirage Full Face

## ResMed Hospital Nasal Mask

ResMed's hospital nasal mask is designed for single patient use. It is available in two sizes—**medium and large**—providing a comfortable fit for a wide range of patients. The masks can be easily discarded, eliminating cross-contamination concerns and costs associated with mask disinfection and sterilization.



ResMed flow generator setting: Standard

**Saving you money...  
saving you time**



### So easy to use, so quick to fit

- Preassembled and ready to use straight from the package
- Color-coded labeling for easy size identification
- Fitting Guide on packaging ensures correct size selection
- Robust one-piece design—no assembly required
- Simple and thoughtful design with minimal adjustment points

### Comfortable and effective therapy

- Effective CO<sub>2</sub> washout
- Reliable therapy—custom designed for the clinical environment
- ResMed sealing technology provides optimal comfort, leak-free seal and reduction of pressure sores
- A broad range of cushion sizes offers a comfortable fit